

System Specifications

For Instrument Manager™ v8.16.x

Each customer situation is unique and system requirements should be discussed with Data Innovations® Sales and Technical teams in advance of any purchase. For budgetary and other planning considerations, Data Innovations has outlined three categories of system specifications to provide guidelines for the type of system that will best suit the intended use of Instrument Manager™ and your production volume. It is important to work with Data Innovations on specific intended use of Instrument Manager and production volumes prior to purchasing computer equipment.

Customer Profile A	Customer Profile B	Customer Profile C
<ul style="list-style-type: none">• Lower volume environments (<1K samples/day)• Instrument Manager used for connectivity only (pass through)• Minimal use of additional IM modules	<ul style="list-style-type: none">• Moderate production volumes (1K-10K samples/day)• Used for connectivity and active lab operations• Data review with Rules for some instruments (other instruments on connectivity only)• Light usage of additional IM modules• Moderate use of Data Archiving	<ul style="list-style-type: none">• Moderate to heavy production volumes. (>10K samples/day)• Heavy use system and lab environment• Significant use of Rules for Autoverification• Multiple users reviewing results in workspaces• Multiple additional IM modules in use• Longer term archiving of information
System Specifications:	System Specifications:	System Specifications:
<ul style="list-style-type: none">• 4 GB of RAM• 80 GB HD (can be multiple drives)• Dual Core Processor• Supported 64-bit OS	<ul style="list-style-type: none">• 8 GB of RAM• 1+ TB HD (can be multiple drives)• Quad Core Processor• Supported 64-bit Server OS	<ul style="list-style-type: none">• 12-16 GB of RAM• 1-5 TB HD (can be multiple drives)• 8 Core Processor• Supported 64-bit Server OS

Report Server System Specifications:

- The Report Server systems specifications, which is for Laboratory Intelligence, should match those of the Primary Production system plus an additional 20% storage for indexing.

Supported Operating Systems

- Windows 7 Professional (64-bit)
- Windows 8.1 Professional (64-bit)
- Windows 10 Professional 64-bit
- Windows Server 2008 R2, Standard Edition (64-bit)
- Windows Server 2012 R2 (64-bit)
- Windows Server 2016

Mirror Set Specifications:

- Disaster Recovery and High Availability both use the mirroring technology. The system specifications should match the system chosen for the Primary Production system.
- Ports 2188 and 443 are needed for Mirroring
- High Availability systems use an Arbiter which have the following specifications:
 - 4 GB of RAM
 - 21 MB of disk space
 - Dual Core Processor
 - Supported 64-bit OS
 - Does not need a dedicated system
 - Must not be installed/running on the active primary or mirror system

End User Specifications

For Instrument Manager™ v8.16.x

Thin-Client Deployment

Supported Operating Systems ^{1,2}

- Same as Primary system with the inclusion of Windows 7 Professional 32-bit. See reverse side

Memory Specifications ³

- 256 MB (16 MB additional memory required on the Base System)

Hard Disk Specifications ³

- Approximately 50 MB for the installation
- 7200 RPM or higher

Other (Hardware) ⁶

- Monitor with minimum 1024 x 768 resolution (1920 x 1080 for Laboratory Intelligence), Keyboard, Mouse, CD/DVD drive for install disk, and an optional bar code scanner

Requirements To Install

- Administrator privileges - Local Administrator logon or local user account access with administrator privileges

Requirements To Use

- Read, write and delete access to the shared Instrument Manager folder on the base Instrument Manager workstation
- Local user access – on the Thin Client Computer workstation

Network Specifications ^{4,5}

- Network interface card installed (NIC)
- TCP/IP networking installed
- Dynamic Host Configuration Protocol (DHCP) or Static IP address
- Port 1972, 57772, and 443 are needed to communicate with Caché
- Remote Procedure Call (RPC services) must be running to configure IM drivers

Supported Browsers

- Microsoft™ Edge
- Microsoft™ Internet Explorer 11

Note: Data Innovations validates using Internet Explorer and Edge, although customers may use other browsers but are responsible for validating the use of the browser.

Browser-Based Deployment (Microsoft Remote Desktop)

- Client computers must have a version of Remote Desktop Connection (RDC) that supports at least Remote Desktop Protocol (RDP) 6.1.
- Minimum 256 MB of additional memory allocated on the server per client PC connected
- ActiveX Client Control must be installed and active.

Notes

¹ Microsoft server class operating systems support a maximum of five (5) Thin Clients out of the box. Additional Thin Clients can be supported by purchasing additional CAL(s) from your authorized Microsoft distributor.

² Microsoft Windows 7 Professional supports a maximum of twenty (20) Thin Clients.

³ The specifications listed are necessary for Instrument Manager Thin Client use. Remember to include operating system and other application software and hardware requirements to ensure full functionality.

⁴ It is recommended that the Primary Instrument Manager system be installed with a Static IP address.

⁵ The InterSystems Caché ODBC driver must be loaded and a System DSN connection must be established if the ODBC Database Access module has been licensed and will be used on the Thin Client workstation. Must use InterSystems' Caché ODBC driver with ODBC Database Access (IM-ODBC-01).

⁶ The Bar Code Scanner to be used with Specimen Storage and Retrieval (SSR), Manual Result Entry (MRE) and SM Workspace within Instrument Manager.

⁷ Other browsers such as Google Chrome, Mozilla Firefox, etc. may be compatible with Microsoft's Remote Desktop Web (RD Web) but have not been tested and their capabilities and/or restrictions are unknown.