



Data Mining Solution

Today's clinical labs are experiencing record growth, which often makes it easy to overlook information that could be critical to ensuring the efficient operation of your laboratory. Since it doesn't look as if the numbers of doctors and patients requiring your services is going to diminish any time soon, it's important to find a way to access key information providing better insight to your lab's operations. No matter how many sites, instruments or customers you serve, Data Innovations' middleware conveniently consolidates all the data in one location for easy retrieval for analysis and reporting.

Imagine the workflow speed with all your data in one location, conveniently consolidated for easy retrieval, analysis and reporting.

Today's clinical labs are experiencing record growth, which often makes it easy to overlook information that could be critical to ensuring the efficient operation of your laboratory. Since it doesn't look as if the numbers of doctors and patients requiring your services is going to diminish any time soon, it's important to find a way to access key information providing better insight to your lab's operations. No matter how many sites, instruments or customers you serve, Data Innovations' middleware conveniently consolidates all the data in one location for easy retrieval for analysis and reporting.

Features of the Data Mining include:

- Ad hoc report generation with any ODBC compliant application such as Microsoft[®] Excel or Crystal Reports[®]
- Automated generation of information and delivery via email or file server on a daily, weekly, monthly basis
- Logical names for database fields for easier data mining
- Ease mundane tasks such as determine autoverification rates, instrument utilization, correlations for instrument or lot comparisons, TAT calculations, and more
- Extract information for EP Evaluator[®], Analyze-it[®] or 'home-grown' QA spreadsheets

Not only do we provide swift and easy access to critical information, at Data Innovations' our unique, powerful Notifier application ties your event notifications—and the recipient—directly to the rules you choose. For example, if you are looking for a sample that meets certain criteria for a study, Notifier can alert you to when that sample is processed.