



# Storage and Retrieval Solution

The workload of the clinical lab has doubled, sometimes tripled, in just a matter of a few years. Data Innovations has solutions for improving efficiency and accuracy in every area, including the creation of a reliable storage and retrieval system, which in this time-sensitive era is more important than ever before. No matter how complicated the footprint of your lab, we have ways to improve workflow efficiency and help you find stored samples with speed and accuracy.

Let us IMpower® your lab to be:

- Vendor Independent
- Work Area Independent
- Sample type independent

## **Other features of our Storage & Retrieval Solution include:**

- Single, networkable sample storage and retrieval solution for your entire enterprise
- Find any specimen in refrigerated storage in under two minutes
- Storage, organization and retrieval of any sample type
- Notify a technician of pending tests or tests being held for verification before storing a sample
- Search by patient ID, sample ID, or patient name with wildcard searches
- Distinguishes different sample types in non-unique sample ID environment (e.g. Same accession number for Chemistry vs. Coagulation)
- Can be integrated with automated specimen storage devices

Time is of the essence in the clinical laboratory; don't waste it. For years to come, you'll be finding specimens in seconds with a system that paid for itself in weeks.