AWS Connectivity Check

Please note the hands-on exercises completed during training are crucial to the learner's success. Inability to access the Virtual Lab will limit learners to the auditory and visual segments of training only.

Test Your Connection to AWS WorkSpaces

- 1. Go to: https://us-east-1.webclient.amazonworkspaces.com/login
- 2. Click Change registration code, enter SLiad+5NUV3B, click Register
- Login with the following credentials: Username: TestConnect Password: Connect1

A successful connection will look like this. No further action is needed.



If you do **not** see the referenced image this means you were not able to connect and will need to contact your IT department to have ports **4195 TCP** and **4195 UDP** opened. These ports are essential to use AWS WorkSpaces. You may also need to adjust your firewall.

Details with additional port requirements can be found by clicking <u>HERE</u>.

Please contact <u>northamerica-traning@datainnovations.com</u> if any questions.