



# Autoverification from Instrument Manager

Time and confidence to focus on what matters most

Today's lab faces greater demands than ever. When you combine a shrinking pool of qualified personnel with increasing workloads and demands for decreased turnaround time, it's more important than ever that labs have the right solutions to help them maximize expertise, streamline operations, and put patient care first.

## Why Autoverification?

- 1 One solution supports all instrument vendors and healthcare platforms, and multiple disciplines
- 2 Autoverification of up to 99% of your test results
- 3 Flexible rule and workflow creation to support unique lab operational needs

Imagine adding an army of virtual lab techs dedicated to alleviating manual review tasks. With Instrument Manager™ Autoverification you can. Autoverification delivers an endlessly flexible rules-writing solution that helps you automate the essential but tedious and time-consuming work of validating every test result, freeing your lab technologists to put their expertise to work by focusing on true exceptions.

## Your lab. Your rules

Lab workflows aren't one-size-fits-all. Lab operations – and the workflows that enable them – are designed to meet the unique requirements of the hospitals, medical professionals, and technology ecosystem that they serve. Autoverification from Instrument Manager is designed by laboratory experts for lab experts.

---

“Instrument Manager Autoverification was a lifesaver for us. We improved workflow efficiency with limited staffing, and improved test turnaround times.”

- Lydia Archuleta

*Clinical Lab Scientist & Middleware Project Manager, Valley Children's Hospital*

## Autoverification from Instrument Manager delivers:

- **Coverage for your total lab** Encompass all of your lab operations with a single solution. Instrument Manager supports a completely open and expansive lab ecosystem, able to integrate with virtually any instrument in your lab – regardless of vendor or discipline. And, it works seamlessly with any LIS, EMR or EHR solution.
- **Your rules, your way** Quickly design and implement customized rules to guide auto-verification of test results across your lab. There's no need to modify your workflows or processes - Autoverification is flexible to adapt to your protocols. And, your lab team can make the rule additions or changes you need, when you need them.
- **Faster results for the people who need them** Anything that bottlenecks lab operations can affect care, delaying results to practitioners and their patients. With Autoverification from Instrument Manager, auto-verified results seamlessly pass through to the hospital's platform, massively speeding delivery of verified results so the right treatment can be applied sooner.
- **Save time and money** Budgets are tight. Time is tighter. Your lab staff are overburdened, but still expected to process more tests, faster, and often without additional help. Autoverification from Instrument Manager eliminates manual work – increasing lab capacity to take on more tests, or to simply complete your existing workload more efficiently. By enabling auto-verification of up to 99% of test results, it's like adding virtual staff without the cost.

Amplify your resources and streamline your workflows with Instrument Manager Autoverification today.

Contact Us

## Key Benefits



Flexible criteria and rule creation to evaluate and take action on every data element in the system



Ability to create custom patient, specimen, and test-level data elements



Customizable results review screens support graphics, previous results, filtering, and color coding



Multi-channel event notifications - including workspace, email, text, and lightpole - for results requiring review



Test environment lets you simulate instrument and/or LIS data to verify rules logic