



Daily & Peer QC Integration with Instrument Manager

Simplify and optimize QC with 3rd party QC program integration

Topics

- I. Understanding Daily and Peer QC Integration from Instrument Manager

Integration with QC program industry leaders

ThermoFisher
S C I E N T I F I C
Thermo Scientific LabLink xL



T E C H N O P A T H
Technopath IAMQC



Bio-Rad QC OnCall™
and Unity Real Time®



CODASY
Codasy QCML*

* Excludes North America

Inadequate QC can have major consequences



Patient
misdiagnosis



Delays in patient
treatment



Inappropriate
treatment



Increased costs due to
retests or unnecessary
further investigations

Cost of avoidable
re-tests


\$200MM USD per year

What if you could reduce QC risk?

- Protect patients and ensure result accuracy
- Meet your lab's tailored workflow needs
- Enhance and accelerate Peer comparison insights



Daily & Peer QC Integration with Instrument Manager



Maximize the value of 3rd party lab QC programs with integrated result delivery and optimized lab workflows


Powerful bi-directional integration



Prevent release of results that fail QC with Autoverification



Deliver QC results to QC program automatically




Receive flags from QC program that can be used to automatically initiate lab workflows

Enhance Peer QC Insights



- Insert QC result into Peer Comparison pool automatically
- Capture every QC result
- Access on-demand comparison with peer labs as QC results are reported

How it works



Daily & Peer QC Requirements

- Specimen Management
- Instrument Manager v8.11
 - Technopath IAMQC Requires IM v8.14
- Instrument Manager Bi-Directional Connection

Data Innovations Services and Support have you covered

92%

of team are medical
technologists

300+

years combined lab
experience

500+

years combined
industry experience

**PLEASE CONTACT DATA INNOVATIONS
SALES**

**northamerica-
sales@datainnovations.com**