



datainnovations.com

Instrument Manager v8.14.10

Overview of the New Features and Functionality

IM v8.14: Specimen Routing & Test Code Mapping

- Specimen Routing (SR)
 - SR Configuration can be edited without stopping SR connection
- IHE Law Support
- Test Code Mapping
 - Improvements when returning ordered test code if multiple occurrences of a test are ordered and have different dilutions or AWOS ID.
- Specimen Management
 - Test Sort in Previous Results Pane
 - Import/Export Color Coding

Demonstrate in the Product

- Specimen Management
 - Run Pane Improvements
 - Sorting Runs
 - Non-Repeatable Data Elements
 - Specimen Warnings Rules and on Release
 - Deleting Specimen

Demonstrate in the Product - *continued*

- Default Color Scheme to improve legibility
- Upgrade Across Hardware Confirmation

Questions?

Thank you for your time!

Thank you for your time!

Questions / Next Steps

Support Slides